

IPG Cantech Process Indicator Tape 143C

EXPOSURE INTERPRETATION:

Steam Sterilization Type 1

Tape is a single use device that is interpreted for appropriate signal color after processing.

Before Cycle: Unexposed 143C




Designed for use in gravity and pre-vacuum sterilizers

Cycle Completed: Exposed 143C



Indicator lines on the tape turn black when exposed to proper steam sterilization conditions.

Lead is not used in the manufacturing of the ink of indicator lines.



Other examples that show acceptable "exposed" results